

Tiani-Spirit - IHE Integration Statement		13 September 2006
Vendor	Product Name	Version
Tiani-Spirit GmbH	Spirit-Archive	1.2
This product implements all of the transactions specified in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options specified below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow (SWF)	Image Manager/Image Archive Performed Procedure Step Manager	none
Patient Information Reconciliation (PIR)	Image Manager/Image Archive Performed Procedure Step Manager	none
Consistent Presentation of Images (CPI)	Image Manager/Image Archive	none
Presentation of Grouped Procedures (PGP)	Image Manager/Image Archive	none
Key Image Note (KIN)	Image Manager/Image Archive	none
Evidence Documents (ED)	Image Manager/Image Archive	none
Simple Image and Numeric Report (SINR)	Report Repository	none
Access to Radiology Information (ARI)	Image Manager/Image Archive Report Repository	none
Audit Trail and Node Authentication (ATNA)	Secure Node Image Manager/Image Archive Performed Procedure Step Manager	encryption
Consistent Time (CT)	Time Client	none
Link to Standard Conformance Statement of the Implementation: www.tiani-spirit.com/documents/is_archive.pdf		
Date of Statement: 13. September 2006		