

Tiani-Spirit - IHE Integration Statement		13 September 2006
Vendor	Product Name	Version
Tiani-Spirit GmbH	Spirit Electronic Health Record	1.0
This product implements all of the transactions specified in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options specified below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Cross-Enterprise Document Sharing (XDS)	Document Registry Document Repository Document Source Document Consumer	none
Cross-Enterprise Sharing of Medical Summaries (XDS-MS)	Document Registry Document Repository Document Source Document Consumer	none
Cross-enterprise Document Sharing for Imaging (XDS-I)	Document Registry Document Repository Document Consumer Imaging Document Source Imaging Document Consumer	none
Patient Identifier Cross-referencing (PIX)	Patient Identity Consumer Patient Identity Cross-reference Manager	none
Patient Demographics Query (PDQ)	Patient Demographics Supplier Patient Demographics Consumer	
Personnel White Pages (PWP)	Personnel White Pages Directory Personnel White Pages Consumer	none
Audit Trail and Node Authentication (ATNA)	Secure Node Audit Record Repository Document Registry Document Repository Document Source Document Consumer Patient Demographics Supplier Patient Demographics Consumer Patient Identity Consumer Patient Identity Cross-reference Manager	encryption
Consistent Time (CT)	Time Client	none
Link to Standard Conformance Statement of the Implementation: www.tiani-spirit.com/documents/is_ehr.pdf		
Date of Statement: 13 September 2006		