

Tiani-Spirit - IHE Integration Statement		13 September 2006
Vendor Tiani-Spirit GmbH	Product Name Spirit Viewer	Version 0.5
This product implements all of the transactions specified in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options specified below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Consistent Presentation of Images (CPI)	Image Display Evidence Creator	none
Evidence Documents (ED)	Image Display Evidence Creator	none
Access to Radiology Information (ARI)	Image Display	none
Audit Trail and Node Authentication (ATNA)	Secure Node Image Display Evidence Creator Document Source Document Consumer Imaging Document Source Imaging Document Consumer	encryption
Consistent Time (CT)	Time Client	none
Cross-Enterprise Document Sharing (XDS)	Document Source Document Consumer	none
Patient Care Coordination (XDS-MS)	Document Source Document Consumer	none
Cross-enterprise Document Sharing for Imaging (XDS-I)	Document Consumer Imaging Document Source Imaging Document Consumer	none
Link to Standard Conformance Statement of the Implementation: www.tiani-spirit.com/documents/is_viewer.pdf		
Date of Statement: 13 September 2006		