

IHE Integration Statement			
Vendor	Product Name	Version	Date
Tiani Sprit "GmbH - Cisco Systems Inc."	EHR_Tiani-Cisco_01		11/30/16

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Access to Radiology Information	Image Display	None
Access to Radiology Information	Image Manager/Archive	None
Audit Trail and Node Authentication	Audit Record Repository	None
Audit Trail and Node Authentication	Secure Application	None
Audit Trail and Node Authentication	Secure Node	None
The Normal Treatment Planning (Simple)	Archive (including RT)	None
The Normal Treatment Planning (Simple)	Dose Displayer	None
Cardiac Catheterization Workflow	Image Display	None
Cardiac Catheterization Workflow	Image Manager/Archive	None
Cardiac Catheterization Workflow	Performed Procedure Step Manager	None
Consistent Presentation of Images	Image Display	None
Consistent Presentation of Images	Image Manager/Archive	None
Consistent Time	Time Client	None
Consistent Time	Time Server	None
Retrieve ECG for Display	Display	None
Echocardiography Workflow	Image Display	None
Echocardiography Workflow	Image Manager/Archive	None
Cardiology Evidence Documents	Image Display	None
Cardiology Evidence Documents	Image Display	Echo Evidence
Cardiology Evidence Documents	Image Manager/Archive	None

Internet address for vendor's IHE information :

Links to Standards Conformance Statements for the Implementation

HL7

DICOM

Link to general information on IHE : www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
Tiani Sprit "GmbH - Cisco Systems Inc."	EHR_Tiani-Cisco_01		11/30/16

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Intravascular Imaging Profile	Image Display	None
Intravascular Imaging Profile	Image Manager/Archive	None
Key Image Note	Image Display	None
Key Image Note	Image Manager/Archive	None
Mammography Image	Image Display	None
Mammography Image	Image Manager/Archive	None
Patient Administration Management	Patient Demographics Consumer	None
Patient Administration Management	Patient Demographics Consumer	Merge
Patient Administration Management	Patient Demographics Consumer	Link / Unlink
Patient Administration Management	Patient Encounter Consumer	Inpatient / Outpatient Encounter Management
Patient Administration Management	Patient Encounter Consumer	None
Patient Demographics Query	Patient Demographics Consumer	None
Patient Demographics Query	Patient Demographics Supplier	None
Patient Demographic Query HL7 V3	Patient Demographics Consumer	None
Patient Demographic Query HL7 V3	Patient Demographics Supplier	None
Patient Identifier Cross-referencing for MPI	Patient Identifier Cross-reference Consumer	None
Patient Identifier Cross-referencing for MPI	Patient Identity Source	None
Patient Identifier Cross-referencing for MPI	Patient Identity Cross-reference Manager	None

Internet address for vendor's IHE information :

Links to Standards Conformance Statements for the Implementation

HL7

DICOM

Link to general information on IHE : www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
Tiani Sprit "GmbH - Cisco Systems Inc."	EHR_Tiani-Cisco_01		11/30/16

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Patient Identifier Cross-Reference HL7 V3	Patient Identifier Cross-reference Consumer	None
Patient Identifier Cross-Reference HL7 V3	Patient Identity Source	None
Patient Identifier Cross-Reference HL7 V3	Patient Identity Cross-reference Manager	None
Radiation Exposure Monitoring	Image Manager/Archive	None
Scheduled Workflow	Image Display	None
Scheduled Workflow	Image Manager/Archive	None
Scheduled Workflow	Performed Procedure Step Manager	None
Cross Community Access	Initiating Gateway	None
Cross Community Access	Responding Gateway	None
Cross-Community Access for Imaging	Imaging Document Consumer	None
Cross-Community Access for Imaging	Imaging Document Source	None
Cross-Community Access for Imaging	Initiating Image Gateway	None
Cross-Community Access for Imaging	Responding Image Gateway	None
Cross-Enterprise Document Sharing for Imaging	Document Consumer	None
Cross-Enterprise Document Sharing for Imaging	Document Registry	None
Cross-Enterprise Document Sharing for Imaging	Document Repository	None

Internet address for vendor's IHE information :	
Links to Standards Conformance Statements for the Implementation	
HL7	
DICOM	
Link to general information on IHE : www.ihe.net	

IHE Integration Statement			
Vendor	Product Name	Version	Date
Tiani Sprit "GmbH - Cisco Systems Inc."	EHR_Tiani-Cisco_01		11/30/16

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Cross-Enterprise Document Sharing for Imaging	Imaging Document Consumer	None
Cross-Enterprise Document Sharing for Imaging	Imaging Document Source	None
Cross-Enterprise Document Sharing for Imaging	Imaging Document Source	Set of DICOM Instances
Cross-Enterprise Document Sharing for Imaging	Imaging Document Source	PDF Report
Cross-Enterprise Document Sharing for Imaging	Imaging Document Source	Text Report
Cross-Enterprise Document Sharing	Document Consumer	None
Cross-Enterprise Document Sharing	Document Registry	None
Cross-Enterprise Document Sharing	Document Repository	None
Cross-Enterprise Document Sharing	Document Source	None

Internet address for vendor's IHE information :

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	

Link to general information on IHE : www.ihe.net